


CyberKnife AQA Report


Left Superior Image



Left Superior Image Coordinates

File Name	C:\RIT\RIT113V4\Demo ...\20160427-LSiris-317.tif
Color Channel	all
Top Orientation	Inferior
Left Orientation	Left
X Centroid Offset (mm)	0.84
Y Centroid Offset (mm)	-0.6
Eccentricity (Beam)	0.19
Eccentricity (Shadow)	0.18

Anterior Superior Image



Anterior Superior Image Coordinates

File Name	C:\RIT\RIT113V4\Demo ...\20160427-ASiris-318.tif
Color Channel	all
Top Orientation	Inferior
Left Orientation	Posterior
X Centroid Offset (mm)	0.51
Y Centroid Offset (mm)	-0.09
Eccentricity (Beam)	0.15
Eccentricity (Shadow)	0.22

Patient Plane Coordinates

X Offset (superior-inferior)	0.35mm
Y Offset (right-left)	-0.84mm
Z Offset (posterior-anterior)	0.51mm
Radial Error	1.04mm
Radial Error Pass/Fail	

Analysis Settings

Baseline Radial Error (mm)	NaN
Radial Error Difference (+/- mm)	1
Ball Size	2
SAD (cm)	80