



CyberKnife Iris Report

Commissioning Image




Commissioning Image Information

| | |
|---------------|---|
| File Name | C:\RIT\RIT113V4\Demo Images\Cyberknife\Iris\20-iris-226.tif |
| Color Channel | all |


Analysis Image**Analysis Image Information**

| | |
|---------------|---|
| File Name | C:\RIT\RIT113V4\Demo Images\Cyberknife\Iris\20fixed263.tif |
| Color Channel | all |

Radial Profiles



Transverse Profiles



Analysis Results

| | |
|--|---------|
| Image Mean FWHM (mm) | 20.126 |
| Commissioning Mean FWHM (mm) | 19.3968 |
| FWHM Difference (mm) | 0.73 |
| Commissioning Radial FWHM (mm) | 19.2989 |
| Commissioning Transverse FWHM (mm) | 19.7956 |
| Image Radial FWHM (mm) | 20.1178 |
| Image Transverse FWHM (mm) | 19.9017 |
| Maximum Radial % Difference (SU) | 9.69 |
| Maximum Transverse % Difference (SU) | 6 |
| FWHM Difference Pass/Fail | Fail |
| Iris Size (mm) | 20 |
| Iris Type (mm) | Iris |
| Maximum FWHM Deviation from Commissioning (mm) | 0.2 |