Assuring MR Image Quality



Has Never Been This Easy

American College of Radiology Accreditation Phantom for Extremity Scanners

Item: ACR-PHE (Small Phantom)

Specifications as designated by the American College of Radiology

The phantom is an MRI accessory used to evaluate the performance of an MRI machine.

Schedule B number 9018.19.0000

The small phantom is a : 4.5" (114.3mm) diameter cylinder, 4.5" (114.3mm) length cylinder, 5.5" (139.70mm) length including level bars.

The various fixtures inside include a: dinosaur, grid, membrane stack, ramp, resolution block, S-N chamber.

This unit may break open if subjected to extremes of temperature. Keep temperature between 60-80°F (15-27°C).



Room 1104, Crawford House, 70 Queen's Road Central, Central, Hong Kong

Phone: +852 2150 1927 Fax: +852 2961 4800 e-mail: cs@jwpacific.com web site: www.jwpacific.com